



Cat.NO. 101081













Store the kit at 15 to 25°C

Genorise Gel DNA Purification Mini Kit PCR Clean Kit 200

For Research Use only

PCR clean-up, Gel extraction

Protocol-at-a-glance

	PCR clean-up	Gel extraction
1 PCR clean-up: Adjust binding condition Gel extraction: Excise DNA fragment / Solubilize gel slice	 100 µL PCR 300µL Buffer GGA 150µL Buffer GGB (Or 300 µL Buffer GP)	  100 mg gel 300 µL Buffer GGA 55°C 5-10 min 150µL Buffer GGB
2 【Optional steps】 Prepare The Column	 	250 µL Buffer GL 10,000 x g 1 min
3 Bind DNA	 	10,000 x g 1 min
4 Wash silica membrane 5 Wash silica membrane	 	500 µL Buffer GW2 700 µL Buffer GW2 10,000 x g 1 min
6 Dry silica membrane		10,000 x g 1 min
7 Elute DNA	 	10–25 µL Buffer GE RT 1 min 10,000 x g 1 min

2.2 Protocol for PCR clean-up

- ① For DNA clean-up from enzymatic reactions using this protocol, two methods to choose from one of the following

@USE Buffer GGA and Buffer GGB

- add 3 volumes of **Buffer GGA**, Add **Buffer GGB** 0.5 volume of **GGA**
For example, if the PCR product is 100µL, add 300µL **Buffer GGA**, Add 150µL **Buffer GGB**
- mix, proceed with **step2.1 ②** (page 6) of the protocol.

@USE Buffer GP

- add 3 volumes of **Buffer GP**
For example, if the PCR product is 100µL, add 300µL **Buffer GP**.
- mix, proceed with **step2.1 ②** (page 6) of the protocol.

3.Troubleshooting

Problem	Possible cause and suggestions
Incomplete Lysis of Agarose Slices	<p><i>Time and temperature</i></p> <ul style="list-style-type: none"> • Check incubation temperature. Depending on the weight of the gel slice, incubation (step ①) can be prolonged up to 20 min. Vortex every 2 min and check integrity of the gel slice. Very large gel slices can be crushed before addition of Buffer GGA.
Low DNA yield	<p><i>Reagents not prepared properly</i></p> <ul style="list-style-type: none"> • Add indicated volume of 95–100% ethanol to Buffer GW2 Concentrate and mix well before use.
	<p><i>Incompletely dissolved gel slice</i></p> <ul style="list-style-type: none"> • Increase time or add another two volumes of Buffer GGA and vortex the tube every 2 minutes during incubation at 50°C. Small pieces of gel are hardly visible and contain DNA that will be lost for purification.
	<p><i>Insufficient drying of the silica membrane</i></p> <ul style="list-style-type: none"> • Centrifuge 5 min at 10,000 x g or incubate column for 2–5 min at 70°C before elution to remove ethanolic Buffer GW2 completely. Ethanolic contaminations are also indicated by gel-loading problems (samples float out of gel slots). remove the spin cup carefully from centrifuge and collection tube and avoid contact of spin cup with flow-through.
	<p><i>Not enough elution buffer</i></p> <ul style="list-style-type: none"> • Especially when larger amounts of DNA (>5µg) are bound, increase elution buffer volume up to 50 µL.
	<p><i>Isolation of large DNA fragments</i></p> <ul style="list-style-type: none"> • Preheat Buffer GE to 70°C, and incubate on the silica membrane at room temperature for 2 min before centrifugation.

② **【Optional steps】 Prepare the DNA Binding Column:**

- Insert one DNA Binding Column into one Collection Tube.
- Pre-equilibrate the silica fleece in the DNA Binding Column by adding 250 μ L **Buffer GL** to the Column before applying sample. Centrifuge for 1 min at 10,000 x g and discard the flowthrough .



250 μ L Buffer GL
10,000 x g 1 min



③ **To bind DNA:**

- Insert the DNA Binding Column into the Collection Tube.
- Transfer the sample from Step① into the Column.
- Insert the DNA Binding Column assembly into a standard tabletop microcentrifuge.
- Centrifuge for at 10,000 x g.
- ▲ The maximum volume of the column reservoir is 750 μ L.
For sample volumes of more than 750 μ L, simply load and spin again.



10,000 x g 1 min



④ **To wash the preparation:**

- Add 500 μ L of **Buffer GW2** to the upper reservoir of the Column.
- Centrifuge for 30 - 60 s at 10,000 x g and discard the flowthrough.



500 μ L Buffer GW2
10,000 x g 1 min



⑤ **To wash the preparation:**

- Combine the DNA binding column with the Collection Tube.
- Add 700 μ L of **Buffer GW2** to the upper reservoir of the Column.
- Centrifuge 1 min at 10,000xg and discard the flowthrough.



700 μ L Buffer GW2
10,000 x g 1 min



⑥ **After discarding the flowthrough liquid:**

- Centrifuge the DNA binding column for additional 1 min at full speed.
- Discard the Collection Tube.
- ▲ The extra centrifugation time ensures removal of residual Wash Buffer.



full speed
1 min



⑦ **To elute the DNA**

- Insert the DNA binding column into a clean, sterile 1.5 ml Microcentrifuge tube.
- Add 10-25 μ L prewarmed Buffer GE to the upper reservoir of the DNA binding column.
- Centrifuge the tube assembly for 1 min at 10,000xg 1min.
- ▲ Elution with 10 μ L (instead of 25 μ L) increases the final DNA concentration in the eluate , but also decreases the overall DNA yield.
- ▲ Warm up the Elution Buffer to 70 $^{\circ}$ C increases the final DNA concentration in the eluate.



10,000 x g
1 min



Table of contents

1. Introduction	4
2. Protocols	
2.1 Protocol for DNA extraction from agarose gels	5-6
2.2 Protocol for PCR clean-up	7
3. Troubleshooting	7

Kit contents

Genorise Gel Extraction Kit PCR Clean Kit	Cat.NO.101081
DNA binding columns	200
Collection Tubes	200
1.5 ml microfuge tubes	200
Buffer GL 【Optional steps】	55ml
Buffer GP (PCR Clean)	100ml
Buffer GGA (Gel Extraction)	100ml
Buffer GGB (Gel Extraction)	60ml
Buffer GW2	72ml+178ml ethanol
Buffer GE	20ml

All buffer store at RT.

1. Introduction

The Purification system, designed for rapid DNA cleanup, includes:

- For direct purification of double- or single-stranded PCR products (100 bp–10 kb) from amplification reactions and DNA cleanup from other enzymatic reactions.
- For extraction of DNA fragments (70 bp – 10 kb) from standard, or low-melt agarose gels in TAE (Tris•acetate/EDTA) or TBE (Tris•borate/ EDTA) buffer and DNA cleanup from enzymatic reactions.

2. Protocol

This protocol is designed to extract and purify DNA of 70 bp to 10 kb from standard or low-melt agarose gels in TAE or TBE buffer. Up to 400 mg agarose can be processed per spin column. This kit can also be used for DNA clean-up from enzymatic reactions. For DNA clean-up from enzymatic reactions using this protocol, add 3 volumes of **Buffer GP** and 1 volume of **isopropanol** to the reaction, mix, and proceed with step 6 of the protocol. Alternatively, use the Min Elute Reaction Clean-up Kit.

Important points before starting

- Add ethanol (95-100%) to Buffer GW2 before use (see bottle label for volume).

2.1 Protocol for DNA extraction from agarose gels

① **I . Excise the DNA fragment from the agarose gel with a clean, sharp scalpel.**

- Minimize the size of the gel slice by removing extra agarose.



II . Weigh the gel slice in a colorless tube. Add 3 volumes of Buffer GGA to 1 volume of gel (100 mg ~ 100µL)

- For example, add 300 µL of **Buffer GGA** to each 100 mg of gel.
- For >2% agarose gels, add 6 volumes of **Buffer GGA**.
- The maximum amount of gel slice per column is 400 mg; for gel slices >400 mg use more than one column.



300µL Buffer GGA/100 mg Gel

III. Incubate at 55 -65°C for 10 min (or until the gel slice has Completely dissolved). To help dissolve gel, mix by vortexing the tube every 2-3 min during the incubation.

- ▲ Solubilize agarose completely. For >2% gels, increase incubation time.

IV. Add Buffer GGB 0.5 volume of GGA to the sample and mix.

- For example, if the agarose gel slice is 100 mg, add 300 µL **Buffer GGA**. Then add 150µL **Buffer GGB**.
- **Buffer GGB** volume= 0.5 volume **Buffer GGA**.